

## Highland County Health Department Laboratory Tests and Fees

Laboratory Test	CPT CODE	FASTING	TUBE	SPIN	SPECIAL INSTRUCTIONS	FEE
<b>Amylase</b>	82150	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>ALT (SGPT)</b>	84460	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>AST (SGOT)</b>	84450	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Basic Metabolic Panel*</b>	80048	Y	Red-SST	Y	Fasting for 12 hours, refrigerate after spin	\$15
<b>Complete Metabolic Panel*</b>	80053	Y	Red-SST	Y	Fasting for 12 hours, refrigerate after spin	\$20
<b>CBC*</b>	85067	N	Lavendar	N	3ml whole blood, invert mix gently, room temp	\$15
<b>CBC w/ Diff</b>	85025	N	Lavendar	N	3ml whole blood, invert mix gently, room temp	\$25
<b>Cholesterol+</b>	82465	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Creatinine-serum</b>	82565	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Fasting blood glucose</b>	82962	Y	Red-SST	Y	Spin within 45 minutes, room temp	\$20
<b>Free T4</b>	84439	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Free T3</b>	84481	N	Red-SST	Y	Separate tube, spin within 45 minutes, Room temp	\$35
<b>Hemoglobin*</b>	85018	N	Lavendar	N	3ml whole blood, invert mix gently, room temp	\$5
<b>Hematocrit*</b>	85014	N	Lavendar	N	3ml whole blood, invert mix gently, room temp	\$5
<b>Hemoglobin A1C*+</b>	83036	Y	Lavendar	N	3ml whole blood, invert mix gently, room temp	\$20
<b>Hepatic Function Panel*+</b>	80076	N	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$15
<b>Hep B sur. Ab</b>	86706	N	Red-SST	Y	Separate Tube, Spin within 45 minutes, REFRIGERATE	\$35
<b>Iron-TIBC</b>	83550	Y	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Lipase</b>	83690	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>Lipid</b>	80061	Y	Red-SST	Y	Spin within 45 minutes, room temp	\$30
<b>Magnesium</b>	83735	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25

<b>Measles</b>	86765	N	Red-SST	Y	Separate tube, Spin within 45 minutes, room temp	\$80
<b>Mumps IgG/IgM</b>	86735					
<b>Rubella IgG</b>	86762					
<b>Potassium</b>	84132	N	Red-SST	Y	Spin within 45 minutes, room temp	\$25
<b>PSA*+</b>	84154	N	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$25
<b>PT with INR*</b>	85610	N	Blue	N	room temp	\$15
<b>Activated PTT with INR</b>	85730	N	Blue	N	room temp	\$25
<b>Renal Function Panel</b>	80069	Y	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$20
<b>Testosterone*</b>	84403	N	Red-SST	Y	Spin within 45 minutes, room temp	\$30
<b>TSH*</b>	84443	N	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$20
<b>Varicella IgG Titer</b>	86787	N	Red-SST	Y	Separate tube, Spin within 45 minutes, room temp	\$35
<b>Vitamin D 25</b>	82306	N	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$45
<b>Vitamin D 1.25</b>	82652	N	Red-SST	Y	Spin within 45 minutes, REFRIGERATE	\$60
<b>Glucose*</b>						\$5
<b>Pregnancy Test*</b>						\$14
<b>Head check</b>						\$0